

High Value...Low Price
High Value...Low Price
High Value...Low Price



GENESIS™

R100 Series

The ROMET® GENESIS R100 Features:

- ***Peace of Mind 3 Year Risk Free Warranty***
- ***No Tools Required***
- ***Simple to Use - no Volume/Pitch to Adjust***
- ***Durable, Lightweight DuPont Delrin Material***
- ***Oral Adapter & Five Clear Straws***
- ***Removable Lanyard***
- ***One 9V Alkaline Battery****
- ***Great for Primary or Back-Up Use***

****See reverse side for additional options***

DISTRIBUTED BY:



SANITAS - BULGARIA Ltd.



A Veteran Owned SBE

25 Tzar Osvoboditel blvd,
1504 Sofia, Bulgaria,
Tel: +359 2 943-38-38,
Tel/Fax: +359 2 944-16-50
e-mail: sales@speech-aids.com,
www.speech-aids.com



GENESIS R100 SERIES OPERATING INSTRUCTIONS

GENESIS R100 Series Operating Instructions

1. Push the *On/Off Button*, you should immediately get a buzz sound – that means the battery has power available for speaking.
2. Place the Romet Electronic Larynx to your neck and practice finding your “sweet spot” to achieve your best vocal tone. **NOTE:** this spot can be many places on either your neck or the side of your mouth; each person’s spot is different depending upon your skin texture and degree of radiation received.

5. Practice – look into a mirror and repeat simple sentences to achieve vocal clarity. Record your voice, play it back to hear how you sound.

Practice.

Do NOT Attempt to twist the barrel or the top head.

GENESIS R100 Series Battery Instructions:

Remove the battery end cap from the unit and insert the 9V battery into the barrel. The polarity of the battery is interchangeable; the + or – battery poles are reversible. Electrical contact will be made automatically when the barrel end cap is screwed onto the end of the barrel.



**Push On
Release Off**



**Hole for Lanyard
Attachment**

Model R100 Series Specials:

R100 (Includes 1 Alkaline Battery, Oral Adapter w/5 Straws, & Lanyard)
R110 (Includes 1 Lithium Battery & Charger, Oral Adapter w/5 Straws, & Lanyard)
R120 (Includes 2 Lithium Batteries & Charger, Oral Adapter w/5 Straws, & Lanyard)

DISTRIBUTED BY:



SANITAS - BULGARIA Ltd.



A Veteran Owned SBE

25 Tzar Osvoboditel blvd,
1504 Sofia, Bulgaria,
Tel: +359 2 943-38-38,
Tel/Fax: +359 2 944-16-50
e-mail: sales@speech-aids.com,
www.speech-aids.com